

art. ERAS

Kit erogatore bidet per vasi AP118GE - AP117G - AP117RG - 5051/WCGE - LK117RGE - AP115SE - AP114RSE - AP113RSE

Kit bidet jet for Wc AP118GE - AP117G - AP117RG - 5051/WCGE - LK117RGE - AP115SE - AP114RSE - AP113RSE

➔ **FLAMINIA.**



dotazione kit supplied



svitare ugello di erogazione
unscrew dispensing nozzle



munirsi di:
_ chiave a brugola 4mm
_ nastro Teflon

equip yourself with:
_ Allen key - 4mm hexagonal
_ Teflon tape



check 1 - foro uscita
check 1 - exit hole



check 2 - foro entrata
check 2 - entry hole



inserire l'erogatore nel canale ceramico
insert the dispenser into the ceramic canal



preparare per centraggio delle componenti
prepare for centering



centrare attraverso la brugola e avvitare
center through the Allen key and screw



serrare con cura
tighten carefully

